

# Data sheet: Freezing Medium for Human ES/iPS Cells

Cat#: RCHEFM001, Lot: R3M121220

For Research Use Only

## ► Properties and biological tests

Appearance	Transparent orange liquid
Aseptic test <sup>(1)</sup>	Negative (day 14)
Mycoplasma <sup>(2)</sup>	Negative
pH <sup>(3)</sup>	7.74
Endotoxin <sup>(4)</sup>	ND

<sup>(1)</sup>Cultured using thioglycolic acid medium for aseptic test and SCD medium

<sup>(2)</sup>Cultured using Myco Trim TC

<sup>(3)</sup>Measured at room temperature

<sup>(4)</sup>Measured using the Kinetic Chromogenic methodology

**ReproCELL Inc.**

<http://www.reprocell.com>

E-mail: [info\\_repro@reprocell.com](mailto:info_repro@reprocell.com)