

Data sheet: Freezing Medium for Human ES/iPS Cells

Cat#: RCHEFM001, Lot: R3M130226

For Research Use Only

► Properties and biological tests

Appearance	Transparent orange liquid
Aseptic test ⁽¹⁾	Negative (day 14)
Mycoplasma ⁽²⁾	Negative
pH ⁽³⁾	7.64
Endotoxin ⁽⁴⁾	0.105 EU/ml

⁽¹⁾Cultured using thioglycolic acid medium for aseptic test and SCD medium

⁽²⁾ Tested by PCR for the detection of mycoplasma.

⁽³⁾ Measured at room temperature

⁽⁴⁾ Measured using the Kinetic Chromogenic methodology

ReproCELL Inc.

<http://www.reprocell.com>

E-mail: info_repro@reprocell.com