

# Data sheet: Freezing Medium for Human ES/iPS Cells

Cat#: RCHEFM001, Lot: R3M131125

For Research Use Only

MFG date 25, Nov, 2013

Expiry date Nov, 2016 \*Guaranteed when stored at -80 deg.C. (Date by use)

## ► Properties and biological tests

Appearance	Transparent orange liquid
Aseptic test <sup>(1)</sup>	Negative (day 14)
Mycoplasma <sup>(2)</sup>	Negative
pH <sup>(3)</sup>	7.79
Endotoxin <sup>(4)</sup>	Not detection

(1) Cultured using thioglycolic acid medium for aseptic test and SCD medium

(2) Tested by PCR for the detection of *Mycoplasma*

(3) Measured at room temperature

(4) Measured using the Kinetic Chromogenic methodology

## ReproCELL Inc.

<http://www.reprocell.com>

E-mail:

➤ Japanese(日本語): [info\\_jp@reprocell.com](mailto:info_jp@reprocell.com)

➤ English(英語): [info\\_en@reprocell.com](mailto:info_en@reprocell.com)