

# Data sheet: Laminin-5

Cat#: RCHEOT004, Lot: 20140306

For Research Use Only

MFG date 6, Mar, 2014

Expiry date March, 2016 \*Guaranteed when stored at -80 deg.C. (Date by use)

## ► Properties

Appearance	Transparent liquid
Concentration	23 µg/mL
Purity (SDS-PAGE)	>95%


## ► Culture test performed using human iPS cells

### Methods

Human iPS cells (201B7; Takahashi, K. et al., *Cell*, 131, 861-72, 2007) were cultured in the ReproFF2 on Laminin-5-coated dish for more than 7 days. Subsequently, the cell morphology was observed and compared with those of the control. Alkaline phosphatase activity was ascertained using the ALP substrate kit III (vector).

### Results

Alkaline phosphatase activity: Qualified



## ReproCELL Inc.

<https://www.reprocell.com>

E-mail:

➤ Japanese (日本語): [info\\_jp@reprocell.com](mailto:info_jp@reprocell.com)

➤ English (英語): [info\\_en@reprocell.com](mailto:info_en@reprocell.com)