

**Cat No. RCHETP002 Lot No. R2M160621****MFG date 21 Jun 2016**

\*Store at -20°C, Performance guaranteed until Exp. Date

**Expiration date Jun 2019**■ **Properties and biological tests**

|                                    |                          |
|------------------------------------|--------------------------|
| <b>Appearance</b>                  | Transparent brown liquid |
| <b>Aseptic test <sup>(1)</sup></b> | Negative (Day 14)        |
| <b>Mycoplasma <sup>(2)</sup></b>   | Negative                 |

<sup>(1)</sup>Cultured using thioglycolic acid medium for aseptic test and SCD medium.

<sup>(2)</sup>Mycoplasma detection performed by PCR analysis.