

Cat No. RCDC101 Lot No. A01QO11**MFG date 11 Oct 2017**

*Store at -20°C, Performance guaranteed until Exp. Date

Expiration date Oct 2020■ **Properties and biological tests**

Appearance	Transparent orange liquid
Mycoplasma ⁽¹⁾	Negative
pH ⁽²⁾	7.34
Osmotic pressure ⁽³⁾	271 mOsm/kg

⁽¹⁾Mycoplasma detection performed by PCR analysis.

⁽²⁾Measured after 12 hours incubation under conditions of 37°C and 5% CO₂.

⁽³⁾Measured using the freezing point depression method.