

**Cat No. RCDH101 Lot No. A01QD08****MFG date 8 Dec 2017**

\*Store at -20°C, Performance guaranteed until Exp. Date

**Expiration date Dec 2020**■ **Properties and biological tests**

<b>Appearance</b>	Transparent orange liquid
<b>Mycoplasma <sup>(1)</sup></b>	Negative
<b>Sterility test <sup>(2)</sup></b>	Pass
<b>pH <sup>(3)</sup></b>	7.64
<b>Osmotic pressure <sup>(4)</sup></b>	359 mOsm/kg

<sup>(1)</sup>Mycoplasma detection performed by PCR analysis.

<sup>(2)</sup>Sterility was confirmed by visual inspection after culturing for a week.

<sup>(3)</sup>Measured after 12 hours incubation under conditions of 37°C and 5% CO<sub>2</sub>.

<sup>(4)</sup>Measured using the freezing point depression method.