

Cat No. RCHEOT002, 003 Lot No. A01RU22

MFG date 22 Jun 2018 *Store at -80°C, Performance guaranteed until Exp. Date
Expiration date Jun 2021

■ Properties and biological tests

Mycoplasma ⁽¹⁾	Negative
Endotoxin	<0.1 EU/mL


⁽¹⁾Mycoplasma detection performed by PCR analysis.

■ Culture test performed using human iPS cells**Methods**

Human iPS cells (201B7; Takahashi, K. et al., *Cell*, 131, 861-72, 2007) were cultured in the Primate ES cell medium containing 5 ng/mL of bFGF in the presence of feeder cells for 4 days. After 4 passages, the cell morphology was observed.
Alkaline phosphatase activity was ascertained using the ALP substrate kit III (VECTOR).

Results

Alkaline phosphatase activity: Qualified



Scale: 250 µm