

Cat No. RCDN001N: 3.0×10^6 cells

Lot No. A01RN07


*Store in liquid nitrogen

■ Properties and biological tests

Viability	93.4 %
-----------	--------

Mycoplasma ⁽¹⁾	Negative
---------------------------	----------

⁽¹⁾Mycoplasma detection performed by PCR analysis.

■ Immunostained images after 14 days culture

Scale: 100 μm

Primary antibodies

Anti-βIII-tubulin (1:500) Covance, MMS-435P
Anti-Tyrosine hydroxylase (TH) (1:500) abcam, ab75875
Anti-Choline acetyltransferase (ChAT) (1:400) Millipore, AB144
Anti-Vesicular glutamine transporter 1 (Vglut1) (1:500) Sigma, V0389
Anti-γ-Aminobutyric acid (GABA) (1:1000) Sigma, A2052