

Data sheet: Freezing Medium for human ES/iPS Cells

Cat#: RCHEFM001, Lot: R3M110118

For Research Use Only

► Properties

Appearance	Transparent orange liquid
Aseptic test ⁽¹⁾	Negative (day 14)
Mycoplasma ⁽²⁾	Negative
pH ⁽³⁾	7.86
Endotoxin ⁽⁴⁾	ND

⁽¹⁾Cultured using thioglycolic acid medium for aseptic test and SCD medium

⁽²⁾Cultured using Myco Trim TC

⁽³⁾Measured at room temperature

⁽⁴⁾Measured using the Kinetic Chromogenic methodology

ReproCELL Inc.

<http://www.reprocell.com>

E-mail: info_repro@reprocell.com