

Data sheet: Freezing Medium for Human ES/iPS Cells

Cat#: RCHEFM001, Lot: R3M131031

For Research Use Only

MFG date 31, Oct, 2013

Expiry date Oct, 2016 *Guaranteed when stored at -80 deg.C. (Date by use)

► Properties and biological tests

Appearance	Transparent orange liquid
Aseptic test ⁽¹⁾	Negative (day 14)
Mycoplasma ⁽²⁾	Negative
pH ⁽³⁾	7.72
Endotoxin ⁽⁴⁾	0.240 EU/mL

(1) Cultured using thioglycolic acid medium for aseptic test and SCD medium

(2) Cultured using Myco Trim TC

(3) Measured at room temperature

(4) Measured using the Kinetic Chromogenic methodology

ReproCELL Inc.

<http://www.reprocell.com>

E-mail:

➤ Japanese (日本語): info_jp@reprocell.com

➤ English (英語): info_en@reprocell.com