

Data sheet: Freezing Medium for Human ES/iPS Cells

Cat#: RCHEFM001, Lot: R3M140122

For Research Use Only

MFG date 22, Jan , 2014

Expiry date Jan, 2017 *Guaranteed when stored at -80 deg.C. (Date by use)

► Properties and biological tests

Appearance	Transparent orange liquid
Aseptic test ⁽¹⁾	Negative (day 14)
Mycoplasma ⁽²⁾	Negative
pH ⁽³⁾	7.45
Endotoxin ⁽⁴⁾	0.21 EU/mL

(1) Cultured using thioglycolic acid medium for aseptic test and SCD medium.

(2) Tested by PCR for the detection of *Mycoplasma*.

(3) Measured at room temperature.

(4) Measured using the Kinetic Chromogenic methodology.

ReproCELL Inc.

<https://www.reprocell.com>

E-mail:

➤ Japanese(日本語): info_jp@reprocell.com

➤ English(英語): info_en@reprocell.com