

**Cat#: RCHEFM003 Lot: B01PU03****MFG date 03 Jul 2016**

\*Store at -20°C, Performance guaranteed until Exp. Date

**Expiration date Jul 2019**■ **Properties and biological tests**

<b>Appearance</b>	Transparent liquid
<b>Aseptic test <sup>(1)</sup></b>	Negative (day 14)
<b>Mycoplasma <sup>(2)</sup></b>	Negative
<b>pH <sup>(3)</sup></b>	7.23
<b>Endotoxin</b>	0.019 EU/mL

<sup>(1)</sup>Cultured using thioglycolic acid medium for aseptic test and SCD medium.

<sup>(2)</sup>Mycoplasma detection performed by PCR analysis.

<sup>(3)</sup>Measured after 12 h incubation under conditions of 37°C and 5% CO<sub>2</sub>.