

Cat#: RCHEOT001 Lot: A01QJ12

For Research Use Only

Expiry date      3 years from the date of manufacture


## ► Culture test performed using human iPS cells

### Methods

Human iPS cells (201B7; Takahashi, K. et al., *Cell*, 131, 861-72, 2007) were cultured in the Primate ES Cell Medium containing 5 ng/mL of bFGF for more than 4 days. Subsequently, the cell morphology was observed and compared with those of the control. Alkaline phosphatase activity was ascertained using the ALP substrate kit III (vector).

### Results

Alkaline phosphatase activity: Qualified



### **ReproCELL Inc.**

<http://www.reprocell.com>

E-mail:

- Japanese (日本語): [info\\_jp@reprocell.com](mailto:info_jp@reprocell.com)
- English (英語): [info\\_en@reprocell.com](mailto:info_en@reprocell.com)