

Cat#: RCHEFM005 Lot: B01RA26**MFG date** 26 Apr 2018

*Store at -20°C, Performance guaranteed until Exp. Date

Expiration date Apr 2021■ **Properties and biological tests**

Appearance	Transparent red liquid
Aseptic test ⁽¹⁾	Negative (day 14)
Mycoplasma ⁽²⁾	Negative
pH ⁽³⁾	7.31
Endotoxin	0.003 EU/mL

⁽¹⁾Cultured using thioglycolic acid medium for aseptic test and SCD medium.

⁽²⁾Mycoplasma detection performed by PCR analysis.

⁽³⁾Measured after 12 h incubation under conditions of 37°C and 5% CO₂.